LightNeb

- IT** > MANUALE ISTRUZIONI D'USO
- EN** > OPERATING INSTRUCTIONS MANUAL
- DE** > BEDIENUNGSANLEITUNG
- FR** > MODE D'EMPLOI
- NL** > VERTALING VAN DE ORIGINELE
GEBRUIKSAANWIJZING
- EL** > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ
- PL** > INSTRUKCJA OBSŁUGI I KONSERWACJI
- ES** > MANUAL DE INSTRUCCIONES DE USO

ITALIANO

ENGLISH

DEUTSCH

FRANÇAIS


NEDERLANDS

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
POLSKI

ESPAÑOL



A

B

C

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Ultrasound aerosol therapy device
VMT electronic aerosol device (Vibrating Mesh Technology).

These user instructions are provided for the device model P0318EM.

INTENDED USE. Medical device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a Doctor.

INDICATIONS FOR USE. Treatment of respiratory diseases. This device is suitable for nebulising drugs in solution (mucolytics in general), drugs in suspension (e.g. corticosteroids, such as beclomethasone dipropionate and budesonide) and saline. Medications must be prescribed by a doctor who has assessed the patient's general condition.

⚠️ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious.

• Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

⚠️ TARGET GROUP OF PATIENTS. Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental, or sensory). It is the medical personnel's responsibility to assess the patient's condition and capabilities in order to determine, when prescribing use of the device, whether the patient is able to operate the aerosol device safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠️ WARNINGS CONCERNING POSSIBLE MALFUNCTIONS. • Should your appliance fail to perform, please contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation.

• Please also refer to the case history of faults and their resolution.

WARNINGS. • Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use. • Always consult your general practitioner for identification of treatment. • Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medicine, dosage and treatment indications. • Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended. • If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual safe for further reference. • Do not expose the device to extreme temperatures. • If the packaging is damaged or opened, contact the distributor or service centre. • The time required to switch from storage to operating conditions is about 2 hours. • The appliance must not be used like a toy. The utmost attention

should be paid when used by children. • In the event of damage to the device after being dropped (e.g. detached parts), do not use the Medical Device but contact an authorised service centre. • Do not use the appliance while driving, or in any other situation where you should not be distracted in order not to endanger yourself, people or pets nearby and the objects around you. • In the event of not using the appliance for a long time, remove the batteries, as they might oxidise and damage the contacts inside the appliance. • Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance to ensure that there is no damage; if it is damaged, immediately take the product to an authorised service centre or your dealer. • The warranty does not cover the batteries or any damage caused by unsuitable, depleted or improperly stored batteries. • It is forbidden to access the device opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user. • Use the medicine as soon as possible once it has been opened and avoid leaving it in the medicine tray; once therapy has ended, do not leave the medicine inside the medicine tray and carry out hygienic preparation operations.

Choking hazard: Some components of the appliance are small enough to be swallowed by children, therefore, keep the appliance out of the reach of children.

Fire hazard: - Do not place the appliance and the batteries near sources of heat, light or in excessively hot and humid premises. Store as described in the chapter on ENVIRONMENTAL CONDITIONS. - THIS device is not suitable for use in the presence of a flammable anaesthetic mixture with oxygen or nitrous oxide. - Never use batteries other than those described in this user manual. - Only use batteries of known makes. Always replace both batteries, do not use old and new batteries together. AA alkaline Ni-MH or Lithium batteries can be used, provided that the latter comply with IEC 60086-4 (if not rechargeable) and IEC 62133 (if rechargeable). The batteries can be recharged with an appropriate charging device, not supplied with this appliance.

Electrocution hazard: - Do not use the device in the shower or bath and do not use it with wet hands. - The appliance does not contain parts that can be repaired by the user. - Do not wash the appliance under running water or by immersion.

Risk of ineffectiveness of the therapy: - Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For further information, please refer to the package leaflet provided by the drug manufacturer. - Use the nebuliser in the correct position, as upright as possible; do not tilt it beyond an angle of 45 degrees, in any direction, to prevent the drug from spilling into the mouth or being excessively dispersed, reducing the effectiveness of the treatment. - Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended. - The average life of the batteries depends on their make. - It is very important to use batteries of the same make and type. - Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used

Risk of infection: - We recommend personal use of the accessories to avoid any risk of infection. - Follow hygienic preparation before each use. - At the end of therapy, do not leave the medicine inside the medicine Tray and carry out the hygienic preparation operations. - If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS. This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

TROUBLESHOOTING		
Problem	Cause	Remedy
The device switches on but does not nebulise or does it poorly.	Malfunction of the medicine tray and mesh unit.	Immerse the medicine tray and mesh unit (2) in a solution of water (50%) and white wine vinegar, rinse and reassemble as described in section 6.3 SANITISATION .
	The medicine does not flow into the nebulising area.	Ensure the small hole in the centre of the medicine tray cover (2B) is not obstructed, if so, clean it.
	The nebulising speed depends on the type of medicine used.	<ul style="list-style-type: none"> • The inhalation treatment time depends on the type of medicine and on the patient's inhaling capacity. • Angle: it must be kept upright when in use. If it is tilted backward, the medicine does not touch the mesh.
	Liquid has seeped out from the mesh	Gently dry the mesh with a cotton swab or shake to remove any droplets from the medicine outlet area.
<ul style="list-style-type: none"> • The blue LEDs flash and the appliance switches off after few seconds. • The appliance switches off. 	<ul style="list-style-type: none"> • No liquid in the tray. • Five minutes have elapsed from switch-on. 	<ul style="list-style-type: none"> • Check for liquid in the medicine tray (2). • Switch the appliance on again.
<ul style="list-style-type: none"> • Light button (1A) turns red and flashes and the appliance does not nebulise or switches off. • The device works intermittently, with blue lights on and does not turn off unless the batteries are removed. 	Batteries are down or depleted.	<ul style="list-style-type: none"> • Replace the batteries (6) if alkaline. • Recharge the Batteries if rechargeable Ni-MH or Lithium.

The device does not switch on.	Power supply problem.	<ul style="list-style-type: none"> • Check the Batteries have been inserted with the correct polarity. • Batteries down or depleted. Replace or recharge.
<p>If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at http://www.flaemnuova.it/it/info/assistenza.</p>		

DISPOSAL

Appliance

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.







Attachment

Accessories such as masks and mouthpiece are to be disposed of as general waste after a disinfection cycle.

Packaging

COMMUNICATION OF SERIOUS EVENTS. Please report any serious incidents in connection with this device to the manufacturer and the competent authority of the Member State in which you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious worsening of a person's state of health.

SYMBOLS ON DEVICE OR PACKAGING		
 Unique device identifier	 Serial number of the device	 Class II device
 Manufacturer	 Type BF applied part	 Direct Current
 Positive polarity	 Negative polarity	 Attention
 Model number	 Phthalate- and bisphenol-free	 Atmospheric pressure limits
 Medical device	 Temperature limits	 Moisture limits
 Production date	 Before use: Caution check instructions for use	
IP21	Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).	
	Medical CE marking ref. regulation 2017/745 EU and subsequent updates	
	Disposal of batteries: The depleted batteries must be disposed of in the appropriate collection containers.	

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES. Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.

TECHNICAL SPECIFICATIONS

Power supply	2.4 - 3Vdc 1.5W 2 x AA alkaline or rechargeable batteries
Power	1.5W
Transducer frequency	108 KHz +/- 3KHz
Noise level (at 1 m)	40 dBA approx.
Appliance dimensions	5.0(L) x 6.2(D) x 13.7(H) cm
Nebuliser weight	(excluding batteries) 115g
Pocket dimensions	19(L) X 15(D) cm
Nebulization speed*	0.30 ml/min approx
Aerosol Output Rate ** (AOR)	0.171 ml/min
Percentage of fill volume emitted in 1 minute **	8.56 %
Aerosol Output ** (AO)	0.697 ml
Particle Size (MMAD)**	3.93 µm
Geometric Standard Deviation (GSD)**	2.30
Respirable fraction < 5 µm **	61.3%
Aerosol fraction < 2 µm **	21.1%

Aerosol fraction from 2 µm to 5 µm **	40.2%
Aerosol fraction > 5 µm **	38.7%
Residual volume **	12.5 µl
Medicine tray capacity	10 ml

(*) The nebulization speed was measured with a 0.9% saline solution at 23°C for weight loss according to Flaem internal procedure. It may vary depending on the nebulization head supplied, the drug and environmental conditions. The delivery speed values may also vary depending on the patient's breathing capacity.

(**) In vitro characterization performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma according to EN ISO 27427 with 4ml, 0.1% Salbutamol in 0.9% NaCl solution. More details are available on request.

APPLIED PARTS

Type BF applied parts are: patient accessories (3, 4, 5)

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature	Between +10°C and +40°C
Relative air humidity	Between 10% and 75%
Atmospheric pressure	Between 69KPa and 106 KPa

Storage and transport conditions:

Ambient temperature	Between -10°C and +35°C
Relative air humidity	Between 10% and 75%
Atmospheric pressure	Between 69KPa and 106 KPa

Note: In order to store the batteries appropriately, when not in use for an extended period of time, extract the batteries from the device and store it in a place in compliance with the above temperatures.

DURATION

The average age of the various components set out below refers to device use for approximately 10 min a day (2 applications of 5 min. each)

1 Appliance body* - 5 years	2 Tank tray and mesh unit* - 1 year
3, 4 and 5 Masks and mouthpiece* - 1 year	* see next page for reference

THE EQUIPMENT INCLUDES (FIG. A, B):

1 - Appliance body Material: Polypropylene	1A – Power button
2- Medicine tray and mesh unit Material: Polypropylene	1B – Battery compartment cover
3- Mouthpiece Material: Silicone	1C – Battery compartment cover hook
4- Paediatric mask Material: Polypropylene + Thermoplastic Elastomers	1D – Battery compartment
5- Adult mask Material: Polypropylene + Thermoplastic Elastomers	2A – Accessory coupling shank
6- 2x1.5V AA Alkaline Batteries	2B – Medicine tray cover
	2C – Tray release buttons
	2D – Tray cover hook

OPERATING INSTRUCTIONS

1 - Inserting the batteries (Fig. D)

1a. Open the battery compartment (**1D**) by pressing the hook of the cover (**1C**) then lift the cover (**1B**) in the direction indicated by the arrows.

1b. Insert the two **AA (6)** batteries, in the **Battery Compartment (1D)**, taking care to adhere to the exact polarity.

1c. Close the **Battery Compartment Cover (1B)** again, and ensure the cover is properly engaged.

2 - Preparation (Fig. E)

ATTENTION!

Prior to every use, the nebuliser and the accessories must be cleaned or disinfected as indicated in sections 5/6. **DISASSEMBLY, HYGIENIC PREPARATION.**

It is recommended to use the nebuliser and its accessories for personal use, to prevent any infections from spreading.

Only use original Flaem accessories.

Accurately wash your hands before preparing the inhalation mixture.

2.1 Lift the **Hook (2D)** of the **Cover (2B)**.

2.2 Completely open the **Medicine tray cover (2B)**.

2.3 Pour the medicine in the tank (max capacity 10ml).

2.4 Close the **Medicine tray cover again (2B)**.

ATTENTION

The medicine must be nebulised after it has been poured into the tank.

3 - Configuration of the appliance (Fig. A)

Apply the desired accessory on the **Accessory Coupling Shank (2A)**:

- **Mouthpiece (3)**, - **Paediatric mask (4)**, - **Adult mask (5)**

If you use the mask accessory, place it on your face as shown in (Fig. F).

4 - Operation (Fig. G)

For effective therapy, after setting up the appliance, sit down in a comfortable and relaxed position and begin. With LightNeb it is possible to choose between two delivery modes:

"ON DEMAND" and **"CONTINUOUS USE"**.

4.1 Press and release button (1A), it turns blue and the appliance starts operating in **"CONTINUOUS USE" mode.**

4.2 Press and release the button (1A) and the appliance switches off.

4.3 Press and hold button (1A) in inhalation and release button (1A) in exhalation to use the **"ON DEMAND" feature.**

4.4 Place the **mouthpiece (3)** onto the mouth. As an alternative, you may use the **paediatric mask (4)** or the **Adult mask (5)**.

4.5 To increase the effectiveness of the therapy, breathe deeply and slowly and after inhaling, hold your breath for a moment, so that the inhaled medicine can settle along the respiratory tract. Then exhale slowly.

4.6 The appliance switches off automatically:

- 4.6a after 5 minutes of operation in “**CONTINUOUS USE**” mode;
- 4.6b instantly when the button is released in “**ON DEMAND**” mode;
- 4.6c if there is no liquid in the tray after the relevant blue LEDs blinking.

4.7 If the appliance is not used for longer than 5 minutes, place it down horizontally as shown in the picture.

4.8 At the end of each use, store the device with all its accessories in a dry and dust-free place.

 Do not operate the appliance without liquid.

5 - Disassembling the appliance (Fig. I)

5.1 Press simultaneously the **unlock buttons (2C)** on the sides of the appliance.

5.2 Extract the medicine tray and mesh unit.

5.3 Clean completely as illustrated in the chapter on **HYGIENIC PREPARATION**.

6 - Hygienic preparation

Before and after each use, the appliance and the accessories must be cleaned and disinfected properly as described below. Otherwise, certain microorganisms might settle in the appliance, leading to an infection hazard. Do not use alcohol or other solvents to clean the device.

6.1 - Cleaning the accessories, medicine tray and mesh:

6.1.1 Washing (Fig. L):

Wash components 2-3-4-5 in warm drinking water with mild dishwashing liquid (non abrasive), then rinse them under the tap with plenty of warm water, to remove any detergent residues.

6.1.2 Drying:

6.1.2.a After washing and rinsing all the components, shake them to remove excess water and fit them back on as described in the following section **REASSEMBLING THE APPLIANCE**.

At this point, to remove the water that settles in the nebulising tray and mesh (2) press button 1A (Fig. G) and wait for the appliance to stop nebulising. After that, switch off the unit.

ATTENTION This operation is very important, because it prevents the deposit of limescale in the micro-holes of the mesh, which might impair proper nebulisation in the subsequent therapeutic application.

6.1.2.b To dry all the components completely, leave them out in a dry place (not in the bathroom, for example) as pictured on the side (Fig. M).

6.2 - Cleaning the Appliance Body

Use a slightly moistened cloth with antibacterial detergent for cleaning the **Appliance Body (1)**.

Then dry it with a soft paper tissue.

ATTENTION

Do not immerse the **Appliance Body (1)** in water and do not wash it under the tap.
Do not use alcohol to clean the device.

6.3- Sanitisation (Fig. H)

6.3.1 Immerse for at least 30' the **Medicine tray and open mesh unit (2)** and accessories 3-4-5 in a solution of 50% water and 50% white wine vinegar. Rinse them with plenty of hot running water.

6.3.2 Drying

See point 6.1.2

6.4-Disinfection

METHODS PARTS TO BE DISINFECTED		method A)	method B)	method C)	method D)
		Milk bottle disinfection appliance with steam (no microwave)	Boiling with demineralised or distilled water	Electrolytic chloroxidiser	Glutaraldehyde based disinfectant (Korsolex basic 4%, 15 min)
2- Medicine tray and mesh		YES weekly max 52 cycles	YES weekly max 52 cycles	NO	YES MAX 300 cycles
3 Mouthpiece		YES MAX 300 cycles	YES MAX 300 cycles	YES MAX 300 cycles	YES MAX 300 cycles
4- 5 Adult and paediatric masks		YES MAX 300 cycles	YES MAX 300 cycles	YES MAX 300 cycles	YES MAX 300 cycles

6.4.1 Method A: Steam disinfection (Fig. N).

Disinfect the accessories listed in Table 1 by using a milk bottle disinfection appliance (**no microwaves**), the treatment should last at least 15 minutes. It is necessary to use demineralised or distilled water to prevent limescale deposits, which might impair the operation of the nebulisation head. Carefully read the operating instructions of the device used with regard to the amount of water required and the method to be followed for disinfection.

6.4.2 Method B: Disinfection by boiling (Fig. O).

Disinfect the accessories listed in the Table 1 by boiling for 10 minutes. It is necessary to use demineralised or distilled water to prevent limescale deposits, which might impair the operation of the nebulisation head.

Before performing a new application, wait for the components to cool and dry.

6.4.3 Method C and D: Chemical disinfection (Fig. P)

6.4.3 (METHOD C) Use of disinfectant with electrolytic chloroxidiser (active principle sodium hypochlorite): immerse accessories 3, 4, 5 in a solution obtained in compliance with the pro-

portions indicated by the liquid manufacturer.

 **ATTENTION** Excluding the medicine tray and mesh unit (2)

6.4.3 (METHOD D) Use of glutaraldehyde type disinfectant: immerse the medicine tray (2) and the accessories (3, 4, 5) in a solution in accordance with the product manufacturer's instructions (e.g. Korsorex basic, 4% for 15 minutes).

6.4.3.a Leave the accessories immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.

6.4.3.b Take the accessories and rinse them with plenty of warm drinking water.

6.4.3.c DRYING: carefully perform all the operations described in section 6.1.2. of the user manual.

6.4.3.d Dispose of the solution according to the disinfectant manufacturer's instructions.

7 - Reassembling the appliance

7.1 Insert the medicine tray and mesh (2) in the appliance case (1) using the appropriate guide (A) (Fig. Q).

Ensure the accessory shank is on the same side as the button.

7.2 Push until the two parts are fully coupled.

7.3 Close the medicine tray cover (**2B**) (Fig. R).

AVAILABLE SPARE PARTS AND ACCESSORIES

 <p>Mouthpiece (in silicone) Model: C007 code: ACO563P</p>	 <p>Paediatric mask (Soft-touch white) Model: B023 code: ACO729P</p>
 <p>Spare part Full medicine tank code: ACO001C</p>	 <p>Adult mask (Soft-touch white) Model: B022 code: ACO730P</p>
<p>Only use original Flaem accessories and spare parts, no liability is accepted if non-original parts or accessories are used.</p>	

Gerät zur Ultraschall-Aerosoltherapie Elektronisches Aerosol VMT-(Vibrating Mesh Technology).

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© 2024 FLAEM NUOVA®
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.