MicroDrop®

Pro2

Inhalationsgerät

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Instruction Manual

Content

	Page
Product description	3
Assembly diagram	4
Aerosol therapy device	6
Warnings and risks	6 – 8
Trouble shooting	9
Disposal	10
Symbols	11
Technical specifications	12 – 13
Duration	14
Operating instructions	15
Speed selector and valve system	16
Hygienic preparation	17
Air filtration	21
Ordering Information	22 – 23

Product description

MicroDrop® Pro2

Professional nebulizer for the therapy of chronic respiratory diseases


Congratulations on your purchase and thank you for choosing our nebulizer system.

Our goal is to fully satisfy customers by offering them high-quality devices for the treatment of respiratory diseases.

The MicroDrop® Pro2 is designed for the treatment of diseases of the lower airways. With an ideal aerosol particle size and short treatment time, MicroDrop® Pro2 helps young and old.

Carefully read these instructions and keep them for future reference.

Visit our website www.mpvmedical.com to see the whole range of **MPV MEDICAL** products and spare parts.

A = # MicroDrop Pro2 plus

B = # RF7-2

Aerosol therapy device

These user instructions are provided for model MicroDrop Pro2 plus and RF7-2 devices. The MicroDrop® Aerosol Therapy device consists of a compressor unit (A), a nebuliser and some accessories (B) and spare parts (D).

INTENDED USE.

Medical device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a doctor.

INDICATIONS FOR USE.

Treatment of respiratory diseases. Medications must be prescribed by a doctor who has assessed the patient's general condition.

CONTRAINDICATIONS.

- The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious.
- Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS.

The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

TARGET GROUP OF PATIENTS.

Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the aerosol safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT.

This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., and at home.

Warnings and risks

WARNINGS CONCERNING POSSIBLE MALFUNCTIONS

- Should your appliance fail to perform, please contact the authorised service centre for clarification.
- The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation.
- Please also refer to the case history of faults and their resolution.

WARNINGS

- Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use.
- Always consult your general practitioner for identification of treatment.
- Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medicine, dosage and treatment indications.
- If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor.
- Keep this manual safe for further reference.
- If the packaging is damaged or opened, contact the distributor or service centre.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near heat sources, in sunlight or in an excessively hot environment.
- The time required to switch from storage to operating conditions is about 2 hours.
- It is forbidden to access the compressor unit opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.

Suffocation risk:

- Some components of the appliance are small enough to be swallowed by children, so keep the appliance out of the reach of children.

Strangulation risk:

Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

Fire risk:

THIS device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.

Risk of electrocution:

- Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance and the power cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer.
- Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable.
- Always keep the power cable away from hot surfaces.
- Never obstruct the ventilation slots located on both sides of the compressor unit.
- Do not handle the compressor unit with wet hands. Do not use the compressor unit in humid environments (e.g. while bathing or showering). Do not immerse the compressor unit in water; if this happens, pull the plug immediately. Do not pull out or touch the compressor unit while immersed in water, unplug it first. Take it immediately to an authorised by the manufacturer service centre or your dealer.

Risk of ineffectiveness of therapy:

- Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For further information, please refer to the package leaflet provided by the drug manufacturer.
- Use the nebuliser in the correct position, as upright as possible; do not tilt the nebuliser beyond an angle of 30 degrees, in any direction, to prevent the drug from spilling into the mouth or being overdispersed, reducing the effectiveness of the treatment.
- Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended.
- Only use the device in a dust-free environment, otherwise therapy may be impaired.
- Do not obstruct or insert objects into the filter and its housing in the unit.
- In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.
- Only use original manufacturer accessories or spare parts, no liability is accepted if non-original parts or accessories are used.

Risk of infection:

- We recommend personal use of the accessories to avoid any risk of infection.
- Follow hygienic preparation before each use. Ensure that tubing and accessories are not stored in the vicinity of other accessories or devices for different therapies (e.g. infusions).
- At the end of therapy, do not leave the medicine inside the nebuliser and proceed with the hygienic preparation operations.
- If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

Risk of injury:

- Do not place the device on a soft supporting surface such as a sofa, a bed or a tablecloth.
- Always operate it on a hard surface clear of any objects.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS

This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit manufacturer website.

CASE HISTORY OF FAULTS AND THEIR RESOLUTION.

Before performing any operation, switch off the appliance and disconnect it from the mains.

Trouble shooting

PROBLEM	CAUSE	REMEDY
The device is not working	The power cable has not been correctly inserted into the power socket	Plug the power cable correctly into the sockets
The device does not nebulise or nebulises poorly	The medicine was not inserted into the nebuliser	Pour the right amount of medicine into the nebuliser
	The nebuliser has not been assembled correctly.	Disassemble and reassemble the nebuliser correctly according to the connection diagram on the cover.
	The nebuliser is obstructed	Hygienic preparation of the nebuliser. Deposits of medicine due to lack of hygienic preparation of the nebuliser impair its efficiency and function. Strictly follow the instructions in the HYGIENIC PREPARATION chapter.
	The accessories are not properly connected to the appliance	Check the correct connection between the air inlet of the appliance and the accessories (see connection diagram on cover).
	The pipe is bent or damaged or twisted	Unwind the pipe and check it for crushing or punctures. Replace it if necessary.
	The air filter is dirty	Replace the filter
The device is louder than usual	The filter is inserted	Insert the filter properly into the housing

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised manufacturer service centre nearest to you.

DISPOSAL

Compressor unit

In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to, separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Nebuliser and accessories

To be disposed of as general waste after a sanitisation cycle.

Packaging

Product box

Product packaging bag.
Tube packaging bag

Heat-shrinkable accessories

Bag packing bag

NOTIFICATION OF SERIOUS EVENTS.

Please report any serious incidents in connection with this device to the manufacturer and the competent authority of the Member State in which you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

SYMBOLS ON DEVICE OR PACKAGING

UDI	Unique device identifier	SN	Serial number of the device
	Class II device		Manufacturer
	Type BF applied part		Production date
LOT	Batch Code		Alternating current
	Attention	PHthalates & BPA FREE	Phthalate- and bisphenol-free
	See instructions for use	#	Model number
MD	Medical device		Temperature limits
	Moisture limits		Atmospheric pressure limits
	Function switch off		
	Function switch on		Before use: Caution check instructions for use
	Distributor	CE 0051	Medical CE marking ref. regulation 2017/745 EU and subsequent updates
IP21	Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).		

Technical specifications

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES

- Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.
- Use the medicine as soon as possible once it has been opened and avoid leaving it in the nebuliser; once therapy has ended, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

COMPRESSOR UNIT SPECIFICATIONS.

Model: MicroDrop Pro2 plus

Power supply:	230V~ 50Hz 210VA
Max. pressure:	3.5 ± 0.5 bar
Air flow to compressor:	14 l/min approx
Noise level (at 1 m):	57 dB (A) approx
Operation:	Continued
Dimension:	20 (L) x 18 (W) x 11 (H) cm

NEBULIZER TECHNICAL SPECIFICATIONS.

Model: RF7-2

RF7 Dual Speed Plus Nebuliser

Minimum drug capacity:	2 ml
Maximum drug capacity:	8 ml

TECHNICAL SPECIFICATIONS OF DEVICE.

Model: MicroDrop Pro2 plus compressor combined RF7-2

Operating pressure:	1.3 bar
(with neb.)	

With Speed selector B2.4

	In pos MAX	in pos Min with valve
Delivery Rate	0,65 ml/min approx ⁽¹⁾	0,29 ml/min approx ⁽¹⁾
<i>(1) The delivery rate was measured with a 0,9% Na Cl saline solution at 23°C by weight loss according to Flaem internal procedure.</i>		
Aerosol Output (AO) (2)	0.417 ml ⁽³⁾	0.429 ml ⁽²⁾
Percentage of fill volume emitted per minute	6.01 % ⁽³⁾	4.64 % ⁽²⁾
Aerosol Output Rate (AOR)	0.120 mg/min ⁽³⁾	0.093 mg/min ⁽²⁾
Particle Size (MMAD)	2.81 µm ⁽³⁾	2.61 µm ⁽²⁾
Geometric Standard Deviation (GSD)	2.29 ⁽³⁾	2.30 ⁽³⁾
Respirable fraction < 5	75.7 % ⁽³⁾	78.2 % ⁽²⁾
Aerosol fraction < 2	34,5 % ⁽³⁾	37,6 % ⁽²⁾
Aerosol fraction from 2 µm to 5	41.2% ⁽³⁾	40.6% ⁽²⁾
Aerosol fraction > 5	24.3% ⁽³⁾	21.8% ⁽²⁾
Residual volume	0.94 ml ⁽³⁾	0.899 ml ⁽²⁾

(2)(3) In vitro characterization according to EN ISO 27427 with 4ml, 0.1% Salbutamol in 0.9% NaCl solution with adult ventilatory pattern. Performance are likely to be different for paediatric or infant populations.

(2) Characterization performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma

*(3) Characterization performed by University of Parma
Further details are available on request.*

APPLIED PARTS

Type BF applied parts are: patient accessories (B3, B4, B5)
Weight: 1.800 kg

ENVIRONMENTAL CONDITIONS

Operating conditions:	Storage and transport conditions:
Ambient temperature: Between +10°C and +40°C	Ambient temperature: Between -25°C and +70°C
Relative air humidity: Between 10% and 95%	Relative air humidity: Between 10% and 95%
Atmospheric pressure: Between 69 KPa and 106 KPa	Atmospheric pressure: Between 69 KPa and 106 KPa

DURATION

Model: MicroDrop Pro2 plus Service life 2000 hours.
(Compressor unit)


Model: RF7-2 The expected average life is 1 year, however it is advisable to replace the nebuliser every 6 months during intensive use (or sooner if the nebuliser is clogged) to ensure maximum therapeutic effectiveness.
(Nebuliser and accessories)

EQUIPMENT AND MATERIAL INFORMATION		
The equipment includes:		Information on materials
A -	Compressor unit - Model: MicroDrop Pro2 plus A1 - Switch, A2 - Air intake, A3 - Air filter A4 - Nebuliser holder, A5 - Carrying handle A6 - Power cable, A7 - Cable compartment porto	
B -	Nebuliser and accessories - Model: RF7-2 B1 - Connecting tube (compressor/nebuliser unit) B2 - RF7 Dual Speed Plus B2.1 - Lower part, B2.2 - Nozzle B2.3 - Upper part, B2.4 - Speed selector with valve	
	B3 - Mouthpiece with valve B3.1 - Expiratory valve	Polypropylene
	B4 - Adult SoftTouch mask	Polypropylene + Thermoplastic Elastomers
	B5 - Paediatric SoftTouch mask	
	B8 - Manual nebulization control	
	IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing the Nebuliser and accessories, perform the same procedure.	
C -	Carrying bag	
D -	Filter	

OPERATING INSTRUCTIONS

Before each use, wash your hands thoroughly and clean your nebulizer as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from drips. This device is suitable for the nebulisation of medicinal substances (solutions and suspensions), for which aerosol treatment is envisaged; such substances must in any case be prescribed by a doctor. In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.

1. Insert the plug of the power cable (A6) into a mains socket corresponding to the voltage of the appliance. It must be positioned in such a way that disconnection from the power grid is not difficult.
2. Insert the nozzle (B2.2) into the upper part (B2.3) by pressing as indicated by the 2 arrows in the 'Connection diagram' in point B2. Insert the Speed selector with valve (B2.4) into the upper part (B2.3) as shown in the 'Connection diagram' in point B2. Pour the medicine prescribed by the doctor into the lower part (B2.1). Close the nebuliser by turning the upper part (B2.3) clockwise.
3. Connect the accessories as shown in the 'Connection Diagram' on the cover.
4. Sit comfortably while holding the nebuliser in your hand, place the mouthpiece to your mouth or use the mask. If you use the mask accessory, place it on your face as shown in the Figure (with or without the use of the elastic band).
5. Switch on the device by pressing the switch (A1) and inhale and exhale deeply; after inhaling, it is recommended to hold your breath for a moment so that the inhaled aerosol droplets can settle. Then exhale slowly.
6. When the application is finished, switch off the device and disconnect it from the mains.

ATTENTION: If, after the therapy session, an obvious deposit of moisture forms inside the tube (B1), detach the tube from the nebuliser and dry it with the pressured air from the compressor; this action prevents possible mould growth inside the tube.


If the connection between compressor and the tube is difficult, please try to simultaneously rotate and insert the end of the tube to the compressor connector and rotate and extract it to remove it.

Speed selector and valve system

Ways of using the 'RF7 Dual Speed Plus' nebuliser with speed selector and valve system. It is professional, fast and suitable for the administration of all types of drugs, including the most expensive ones, even in patients with chronic diseases. Thanks to the geometries of the internal ducts of the RF7 Dual Speed Plus nebuliser, a suitable and active particle size was obtained for treatment down to the lower respiratory tract.

The nebulizer can deliver the drug in 2 ways:

In pos MAX

To speed up inhalation therapy, position the speed selector with valve (B2.4) by pressing the MAX marking with your finger.


In pos Min with valve


To make the inhalation therapy more effective, position the speed selector with valve (B2.4) by pressing with a finger on the side opposite the Max marking, in this case you will obtain an optimal administration of the drug while minimising its dispersion in the surrounding environment, thanks to the valve system with which the ampoule, mouthpiece and mask are equipped.

Use of nebulization manual control

To achieve continuous nebulisation action you should not use the manual nebulisation control (B8), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (B8) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited.

To disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

Hygienic preparation

Switch off the device before each hygienic preparation operation and disconnect it from the mains.

Compressor unit (A) and tube exterior (B1)

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Nebuliser and accessories

Open the nebuliser by turning the upper part (B2.3) anticlockwise, detach the nozzle (B2.2) and the speed selector with valve (B2.4) from the upper part (B2.3) as shown in the "Connection diagram" in B2. Then proceed according to the instructions below.

Sanitisation

Before and after each use, sanitise the nebuliser and accessories by choosing one of the methods provided in the table and described below.

method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

method B: Sanitise the accessories in the dishwasher with the hot cycle (70°C).

method C: Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection

After sanitising the nebuliser and accessories, disinfect them using one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with lukewarm drinking water.
- Dispose of the solution according to the disinfectant manufacturer's instructions.

method B: Disinfect accessories by boiling them in water for 10 minutes; use demineralised or distilled water to avoid hard water deposits.

method C: Disinfect the accessories in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

If you also want to carry out sterilisation, skip to the Sterilisation section. After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Mechanical Cleaning and Disinfection

a. Equipment

- Cleaning and disinfection device complying to DIN EN ISO 15883
- Cleaning agent
- Neutralizing agent (if necessary)

b. Procedure

- Place the nebulizer and the accessories on the cleaning device tray so that they can be cleaned in the most effective way
- Select a suitable cleaning program and agent, ref. to Validation paragraph.

c. Validation

- Miele PG8582 cleaner & disinfectant device, program VARIO TD
- Standard baskets made by Miele (E486 basket recommended)
- Neodisher Mediclean forte 0.5% and Neodisher Z for neutralization (from Dr. Weigert)

Sterilisation

Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen. The sterilisation procedure was validated in accordance with ISO 17665-1.


At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Table of planned methods / patient accessories								
Method	B2.1	B2.2	B2.3	B2.4	B3	B4	B5	B8
HYGIENIC PREPARATION IN THE HOME								
Sanitisation								
method A	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓	✓
Disinfection								
method A	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300
method C	300	300	300	300	300	300	300	300
HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL								
Disinfection								
method A	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300
Sterilisation								
	300 Assembled nebuliser				300	300	300	300
	✓: planned \: not planned 300: MAX 300TIMES							

AIR FILTRATION

The compressor is equipped with an air filter (A3) to be replaced with a replacement filter when it is dirty or changes colour. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance.

The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters.

To replace the filter, pull it out as shown in the figure.

The filter is designed so that it is always fitted in its housing. Do not replace the filter during use. **Only use original manufacturer accessories or spare parts, no liability is accepted if non-original parts or accessories are used.**

Ordering Information

MicroDrop® Pro2

Professional nebulizer for the therapy of chronic respiratory diseases

REIMBURSEMENT NUMBER 14.24.01.0121

Including: Compressor, RF7 plus nebulizer, Soft-children and -adult mask, mouthpiece, tubing, manual nebulisation control, carrying bag

ARTICLE	PZN**	NO.
MicroDrop® Pro2	03204625	M51804-00

MicroDrop® profi nebulizer kit universal

Including: RF7 plus nebulizer, mouthpiece with exhalation valve, Soft-children and -adult mask, tubing, compressor filter, manual

PZN**	NO.
00347465	M53700-10

MicroDrop® Pro2 Year-Set

Including: RF7 plus nebulizer, mouthpiece with exhalation valve, tubing, compressor filter


PZN**	NO.
03071182	M50020-45300

MicroDrop® RF7 plus profi nebulizer

with mouthpiece and exhalation valve

PZN**	NO.
00348507	M50020-43900

MicroDrop® RF6 plus mouthpiece

with exhalation valve

PZN**	NO.
00023001	M50020-43200

MicroDrop® Pro2 nozzle (spare part)

for RF7 plus nebulizers

PZN**	NO.
03954496	M50020-44500

MicroDrop® RF7 plus soft-mask adults

silver / transparent with rubber band

PZN**	NO.
00347956	M50020-44100

MicroDrop® RF7 plus soft-mask children

silver / transparent with rubber band

PZN**	NO.
00348051	M50020-44000

MicroDrop® tubing 2m

with 2 fits-all connectors phthalate and BPA free

PZN**	NO.
09937470	M50020-42400

MicroDrop® 2 Filter compressor (spare part)

1 piece for new compressor series

PZN**	NO.
09937530	M50010-42800

MicroDrop® Pro2 manual nebulisation control (spare part)


PZN**	NO.
03954639	M50010-43400

MicroDrop® 2 carrying bag (spare part)

blue

PZN**	NO.
09937553	M50010-42900

Hersteller / MANUFACTURER:
Designed und entwickelt von:
DESIGNED AND ENGINEERED BY:

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FLAEM NUOVA S.p.A., Via Colli Storici, 221-223-225, 25015 S.
FLAEM NUOVA S.p.A.,